

## Ampoules for steam sterilization



Description	Packaging	Ref.
<b>Steril Control GST E6 ampoules</b> Ampoules with nutrient Medium, pH indicator and spores of <i>G. stearothermophilus</i> (1-5x10 <sup>6</sup> CFU/ampoule).	20 Ampoules 4 mL	91210
	50 Ampoules 4 mL	91050
<b>Steril Control GST E5 ampoules</b> Ampoules with nutrient Medium, pH indicator and spores of <i>G. stearothermophilus</i> (1-5x10 <sup>5</sup> CFU/ampoule).	50 Ampoules 1 mL + 10 negative controls	91040
	20 Ampoules 4 mL	91251
<b>Steril Control GST E4 ampoules</b> Ampoules with nutrient Medium, pH indicator and spores of <i>G. stearothermophilus</i> (1-5x10 <sup>4</sup> CFU/ampoule).	50 Ampoules 4 mL	91052
	20 Ampoules 4 mL	91253
	50 Ampoules 4 mL	91054

