


### Flexible Contact Slide for microbial screening of surfaces (ISO 18593)



Description	Packaging	Ref.
<b>CONTACT SLIDE 1</b> (PCA + TTC + NEUTRALIZING / V.R.B.G. AGAR + NEUTRALIZING) Total bacterial count, enterobacteria, with inactivation of disinfectants.	20 slides	525262
	120 slides	53526
<b>CONTACT SLIDE 2</b> (PCA + TTC + NEUTRALIZING / ROSE BENGAL CAF AGAR + NEUTRALIZING) Total bacterial count, yeasts and moulds, with inactivation of disinfectants.	20 slides	525272
	120 slides	53527
<b>CONTACT SLIDE 3</b> (PCA+TTC + NEUTRALIZING / VOGEL JOHNSON AGAR) Total bacterial count with inactivation of disinfectants, staphylococci.	20 slides	525302
	120 slides	53530
<b>CONTACT SLIDE 4</b> (CETRIMIDE AGAR + NEUTRALIZING / ROSE BENGAL CAF AGAR + NEUTRALIZING) <i>Pseudomonas</i> spp., yeasts and moulds.	20 slides	525282 #
	120 slides	53528 #
<b>CONTACT SLIDE 5</b> (V.R.B.G. AGAR + NEUTRALIZING / BILE ESCULIN AGAR + NEUTRALIZING) Enterobacteria, fecal streptococci, with inactivation of disinfectants.	20 slides	525312
	120 slides	53531
<b>CONTACT SLIDE 8</b> (PCA + TTC + NEUTRALIZING / V.R.B.L. AGAR + NEUTRALIZING) Total bacterial count, coliforms, with inactivation of disinfectants.	20 slides	525622
	120 slides	53562
<b>CONTACT SLIDE 9</b> (PCA + TTC + NEUTRALIZING / PCA + TTC + NEUTRALIZING) Total bacterial count with inactivation of disinfectants.	20 slides	525322
	120 slides	53532
<b>CONTACT SLIDE 10</b> (ROSE BENGAL CAF AGAR + NEUTRALIZING / VOGEL JOHNSON AGAR) Yeasts and moulds with inactivation of disinfectants, staphylococci.	20 slides	525452 #
	120 slides	53545 #
<b>CONTACT SLIDE Chrom 1</b> (Chromatic™ E.coli-coliform / Chromatic™ Salmonella) Chromogenic selective media for detecting <i>E. coli</i> , coliforms and <i>Salmonella</i> spp. on surfaces or in food and waters.	20 slides	525292 #
	120 slides	53529 #
<b>CONTACT SLIDE Chrom 2</b> (Chromatic™ E.coli-coliform / PCA+TTC + NEUTRALIZING) Chromogenic selective medium for detecting <i>E. coli</i> , coliforms and for total bacterial count.	20 slides	525392 #
	120 slides	53539 #
<b>CONTACT SLIDE Chrom 3</b> (Chromatic™ E.coli-coliform / BAIRD PARKER AGAR + NEUTRALIZING) Chromogenic selective medium for detecting <i>E. coli</i> , coliforms and <i>Staphylococcus aureus</i> with inactivation of disinfectants.	20 slides	525462 #
	120 slides	53546 #
<b>CONTACT SLIDE TSA + NEUTRALIZING / TSA + NEUTRALIZING</b> Total bacterial count with inactivation of disinfectants.	20 slides	525472 #
	120 slides	53547 #
<b>CONTACT SLIDE TSA + NEUTRALIZING / ROSE BENGAL CAF AGAR + NEUTRALIZING</b> Total bacterial count with inactivation of disinfectants, yeasts and moulds.	20 slides	525482 #
	120 slides	53548 #