

<b>Media for susceptibility testing of mycobacteria</b>				# Product with minimum request
Description	µg / mL	Packaging	Ref.	
LOWENSTEIN JENSEN + AMIKACIN	5	10 tubes	34127 #	
LOWENSTEIN JENSEN + CAPREOMYCIN	10	10 tubes	34138/1 #	
LOWENSTEIN JENSEN + CYCLOSERINE	30	10 tubes	34137/1 #	
LOWENSTEIN JENSEN + ETHAMBUTOL	2	10 tubes	34126/1 #	
LOWENSTEIN JENSEN + ETHAMBUTOL	4	10 tubes	34126/2 #	
LOWENSTEIN JENSEN + ETHAMBUTOL	5	10 tubes	34126/3 #	
LOWENSTEIN JENSEN + ETHIONAMIDE	10	10 tubes	34132/1 #	
LOWENSTEIN JENSEN + ETHIONAMIDE	20	10 tubes	34132/2 #	
LOWENSTEIN JENSEN + ETHIONAMIDE	30	10 tubes	34132/3 #	
LOWENSTEIN JENSEN + ISONIAZIDE	0.1	10 tubes	34123 #	
LOWENSTEIN JENSEN + ISONIAZIDE	0.2	10 tubes	34123/1 #	
LOWENSTEIN JENSEN + ISONIAZIDE	1	10 tubes	34123/2 #	
LOWENSTEIN JENSEN + ISONIAZIDE	5	10 tubes	34123/3 #	
LOWENSTEIN JENSEN + KANAMYCIN	10	10 tubes	34143/1 #	
LOWENSTEIN JENSEN + KANAMYCIN	20	10 tubes	34143/2 #	
LOWENSTEIN JENSEN + KANAMYCIN	30	10 tubes	34143/3 #	
LOWENSTEIN JENSEN + LEVOFLOXACIN	2	10 tubes	34146/1 #	
LOWENSTEIN JENSEN + NICOTINAMIDE	10	10 tubes	34135/1 #	
LOWENSTEIN JENSEN + NICOTINAMIDE	20	10 tubes	34135/2 #	
LOWENSTEIN JENSEN + NICOTINAMIDE	30	10 tubes	34135/3 #	
LOWENSTEIN JENSEN + OFLOXACIN	5	10 tubes	34128/1 #	
LOWENSTEIN JENSEN + OFLOXACIN	10	10 tubes	34128/2 #	
LOWENSTEIN JENSEN + OFLOXACIN	25	10 tubes	34128/3 #	
LOWENSTEIN JENSEN + PAS	0.5	10 tubes	34129/3 #	
LOWENSTEIN JENSEN + PAS	1	10 tubes	34129/1 #	
LOWENSTEIN JENSEN + PAS	10	10 tubes	34129/2 #	
LOWENSTEIN JENSEN + PEFLOXACIN	2	10 tubes	34136 #	
LOWENSTEIN JENSEN + PIRAZINAMIDE	5	10 tubes	34124/1 #	
LOWENSTEIN JENSEN + PIRAZINAMIDE	15	10 tubes	34124/2 #	
LOWENSTEIN JENSEN + PIRAZINAMIDE	20	10 tubes	34124/3 #	
LOWENSTEIN JENSEN + PIRAZINAMIDE	200	10 tubes	34124/4 #	
LOWENSTEIN JENSEN + PYRUVATE	0.2%	10 tubes	34144 #	
LOWENSTEIN JENSEN + RIFABUTIN	10	10 tubes	34130/1 #	
LOWENSTEIN JENSEN + RIFABUTIN	30	10 tubes	34130/2 #	
LOWENSTEIN JENSEN + RIFABUTIN	50	10 tubes	34130/3 #	
LOWENSTEIN JENSEN + RIFAMPICIN	5	10 tubes	34121/1 #	
LOWENSTEIN JENSEN + RIFAMPICIN	10	10 tubes	34121/2 #	
LOWENSTEIN JENSEN + RIFAMPICIN	15	10 tubes	34121 #	
LOWENSTEIN JENSEN + RIFAMPICIN	20	10 tubes	34121/6 #	
LOWENSTEIN JENSEN + RIFAMPICIN	25	10 tubes	34121/3 #	
LOWENSTEIN JENSEN + RIFAMPICIN	40	10 tubes	34121/5 #	
LOWENSTEIN JENSEN + RIFAMPICIN	50	10 tubes	34121/4 #	
LOWENSTEIN JENSEN + RIFAPENTIN	9	10 tubes	34122 #	
LOWENSTEIN JENSEN + STREPTOMYCIN	4	10 tubes	34125/1 #	
LOWENSTEIN JENSEN + STREPTOMYCIN	10	10 tubes	34125/2 #	
LOWENSTEIN JENSEN + STREPTOMYCIN	25	10 tubes	34125/3 #	
LOWENSTEIN JENSEN + TCH	2	10 tubes	35148 #	

*Different antimicrobial agents and concentrations can be produced on specific request.*