

Poliflex coupons for hydrogen peroxide sterilization

Product with minimum request



Description	Packaging	Ref.
OXI CONTROL E6 PLX Poliflex coupons with spores of <i>G.stearothermophilus</i> (1-5x10 ⁶ CFU/coupon) and tubes of nutrient broth with pH indicator.	20 coupons 20 tubes	91085
OXI CONTROL E5 PLX Poliflex coupons with spores of <i>G.stearothermophilus</i> (1-5x10 ⁵ CFU/coupon) and tubes of nutrient broth with pH indicator.	20 coupons 20 tubes	91086 #