


Ready culture media in 55 mm contact plates

Product with minimum order



Description		Packaging	Ref.
BAIRD PARKER AGAR + NEUTRALIZING Selective medium for <i>Staphylococcus aureus</i> isolation with inactivation of disinfectants.		20 plates	15331
CETRIMIDE AGAR Selective medium for <i>Pseudomonas aeruginosa</i> isolation.	E.P.	20 plates	15332
LEGIONELLA AGAR (B.C.Y.E.) Non-selective medium for <i>Legionella</i> spp. isolation.	ISO 11731-1	20 plates	15334
LEGIONELLA AGAR (G.V.P.C.) Selective medium for <i>Legionella</i> spp. isolation.	ISO 11731-1	20 plates	15376
MAC CONKEY AGAR Selective medium for gram-negatives isolation.	E.P.	20 plates	15329
MANNITOL SALT AGAR Selective medium for staphylococci isolation.	E.P.	20 plates	15328
m-FAECAL COLIFORM AGAR Selective medium for fecal coliforms isolation and enumeration.		20 plates	15368
PLATE COUNT AGAR Medium for surfaces and air total bacterial count.	ISO 4833	20 plates	15325
PLATE COUNT AGAR + NEUTRALIZING Medium for surfaces total bacterial count with inactivation of disinfectants.	ISO 4833	20 plates	15336
PLATE COUNT AGAR + NEUTRALIZING (IRRADIATED) Medium for surfaces total bacterial count with inactivation of disinfectants.	ISO 4833	20 plates	15336S
PLATE COUNT AGAR + TTC + NEUTRALIZING Medium for surfaces total bacterial count with inactivation of disinfectants.	ISO 4833	20 plates	15360
ROSE BENGAL AGAR + CAF Selective medium for yeasts and moulds isolation and enumeration.	APHA	20 plates	15374
R2 AGAR Medium for heterotrophic microorganisms count in drinking water.	E.P.	20 plates	15354 #
SABOURAUD DEXTROSE AGAR Medium for yeasts and moulds isolation.	E.P.	20 plates	15327
SABOURAUD AGAR + NEUTRALIZING (IRRADIATED) Medium for yeasts and moulds isolation with inactivation of disinfectants.		20 plates	15323S
SABOURAUD CAF AGAR (IRRADIATED) Selective medium for yeasts and moulds isolation.		20 plates	15380S
SABOURAUD CAF AGAR + NEUTRALIZING Selective medium for yeasts and moulds isolation with inactivation of disinfectants.		20 plates	15365
SLANETZ BARTLEY AGAR + NEUTRALIZING Selective medium for fecal streptococci enumeration and isolation with inactivation of disinfectants.		20 plates	15362
TRYPTIC SOY AGAR (IRRADIATED) Medium for surfaces and air total bacterial count.	E.P.	20 plates	15340S
TRYPTIC SOY AGAR + NEUTRALIZING (IRRADIATED) Medium for surfaces and air total bacterial count with inactivation of disinfectants.	E.P.	20 plates	15339S
TRYPTIC SOY AGAR + 0.1% PENASE + NEUTRALIZING Medium for surfaces and air total bacterial count with inactivation of disinfectants penicillins.		20 plates	15330 #
VIOLET RED BILE GLUCOSE AGAR Selective medium for enterobacteria isolation and enumeration.	ISO 21528-2 E.P.	20 plates	15375



Ready culture media in 55 mm contact plates

Product with minimum order



Description		Packaging	Ref.
VIOLET RED BILE LACTOSE AGAR Selective medium for coliforms isolation and enumeration.	ISO 4832	20 plates	15326
VIOLET RED BILE GLUCOSE AGAR + NEUTRALIZING Selective medium for <i>Enterobacteriaceae</i> isolation and enumeration with inactivation of disinfectants.	ISO 21528-2	20 plates	15382
VIOLET RED BILE LACTOSE AGAR + NEUTRALIZING Selective medium for coliforms isolation and enumeration with inactivation of disinfectants.	ISO 4832	20 plates	15361