

Ready culture media in three sectors petri dishes

Product with minimum order



Description	Packaging	Ref.
BAIRD PARKER / BIGGY / MAC CONKEY Selective media for <i>Staphylococcus aureus</i> , <i>Candida</i> spp. and gram-negatives isolation.	20 plates 100 plates	13013 # 13013* #
BAIRD PARKER / BILE ESC. / SABOURAUD CAF Selective media for <i>Staphylococcus aureus</i> , group D streptococci and fungi isolation.	20 plates 100 plates	13602 # 13602* #
BAIRD PARKER / MAC CONKEY / SABOURAUD CAF Selective media for <i>Staphylococcus aureus</i> , gram-negatives and fungi isolation.	20 plates 100 plates	13371 # 13371* #
CHOCOLATE / COLUMBIA CNA / THAYER MARTIN Media for fastidious bacteria, streptococci, staphylococci, pneumococci and pathogenic neisseriae isolation.	20 plates 100 plates	13014 # 13014* #
CLED / ENTEROCOCCO / MAC CONKEY Media for total bacterial count in urine and isolation of enterococci and gram-negatives.	20 plates 100 plates	13614 # 13614* #
CLED / MAC CONKEY / CETRIMIDE Media for total bacterial count and for Gram-negative bacteria and <i>Pseudomonas aeruginosa</i> detection.	20 plates 100 plates	13019 # 13019* #
GARD. V. / ROGOSA / THAYER MARTIN Media for <i>Gardnerella vaginalis</i> , lactobacilli, pathogenic neisseriae isolation.	20 plates 100 plates	13345 # 13345* #
MAC CONKEY / V. JOHNSON / SABOURAUD Media for gram-negatives, <i>Staphylococcus aureus</i> and fungi isolation.	20 plates 100 plates	13480 # 13480* #