

Ready culture media in two sectors petri dishes

Product with minimum order



Description	Packaging	Ref.
BAIRD PARKER / MAC CONKEY Selective media for <i>Staphylococcus aureus</i> and gram-negatives isolation.	20 plates 100 plates	18500 # 18500* #
BAIRD PARKER / SABOURAUD CAF Selective media for <i>Staphylococcus aureus</i> , yeasts and moulds isolation.	20 plates 100 plates	18390 # 18390* #
BRILLIANT GREEN / S.S. Selective media for <i>Salmonella</i> spp. isolation.	20 plates 100 plates	18012 # 18012* #
CLED / MAC CONKEY Media for total bacterial count and gram-negatives isolation in urine.	20 plates 100 plates	18502 # 18502* #
COLUMBIA C.N.A. / CHOCOLATE Selective medium for staphylococci, streptococci, pneumococci and non selective medium for fastidious bacteria isolation.	20 plates 100 plates	18507 # 18507* #
COLUMBIA C.N.A. / GARDNERELLA V. Selective media for streptococci, staphylococci, pneumococci and <i>Gardnerella vaginalis</i> isolation.	20 plates 100 plates	18422 # 18422* #
COLUMBIA C.N.A. / MAC CONKEY Selective media for streptococci, staphylococci, pneumococci and gram-negative isolation.	20 plates 100 plates	18327 # 18327* #
D.T.M. / SABOURAUD Selective media for dermatophytes, yeasts, moulds isolation.	20 plates 100 plates	18595 # 18595* #
GARDNERELLA V. / THAYER MARTIN Selective media for <i>Gardnerella vaginalis</i> and pathogenic neisseriae isolation.	20 plates 100 plates	18379 # 18379* #
HEKTOEN ENTERIC / YERSINIA Selective media for enteric pathogens and <i>Yersinia enterocolitica</i> isolation.	20 plates 100 plates	18391 # 18391* #
HEKTOEN ENTERIC / S.S. Selective media for enteric pathogens isolation.	20 plates 100 plates	18503 # 18503* #
MAC CONKEY / SABOURAUD CAF Selective media for coliforms and pathogenic fungi detection.	20 plates	18018 #
MAC CONKEY / S.S. Selective media for enteric pathogens isolation.	20 plates 100 plates	18505 # 18505* #
MAC CONKEY / T.S.A. BLOOD Media for gram-negatives and fastidious bacteria isolation.	20 plates 100 plates	18380 # 18380* #