

Strips for dry heat or ethylene oxide sterilization		# Product with minimum request
Description	Packaging	Ref.
STRIP CONTROL BAT E8 Strips with spores of <i>B. atrophaeus</i> (1-5x10 ⁸ CFU/strip) and tubes of nutrient broth with pH indicator.	20 strips 20 tubes	91059 #
STRIP CONTROL BAT E7 Strips with spores of <i>B. atrophaeus</i> (1-5x10 ⁷ CFU/strip) and tubes of nutrient broth with pH indicator.	20 strips 20 tubes	91062 #
STRIP CONTROL BAT E6 Strips with spores of <i>B. atrophaeus</i> (1-5x10 ⁶ CFU/strip) and tubes of nutrient broth with pH indicator.	20 strips 20 tubes	91063
	100 strips	91067
STRIP CONTROL BAT E5 Strips with spores of <i>B. atrophaeus</i> (1-5x10 ⁵ CFU/strip) and tubes of nutrient broth with pH indicator.	20 strips 20 tubes	91064 #
STRIP CONTROL BAT E4 Strips with spores of <i>B. atrophaeus</i> (1-5x10 ⁴ CFU/strip) and tubes of nutrient broth with pH indicator.	20 strips 20 tubes	91065 #
STRIP CONTROL BAT E3 Strips with spores of <i>B. atrophaeus</i> (1-5x10 ³ CFU/strip) and tubes of nutrient broth with pH indicator.	20 strips 20 tubes	91066 #