



<b>Strips for irradiation sterilization</b>		# Product with minimum request
Description	Packaging	Ref.
<b>STRIP CONTROL PUM E7</b> Strips with spores of <i>B. pumilus</i> (1-5x10 <sup>7</sup> CFU/strip) and tubes of nutrient broth with pH indicator.	20 strips 20 tubes	91071 #
<b>STRIP CONTROL PUM E6</b> Strips with spores of <i>B. pumilus</i> (1-5x10 <sup>6</sup> CFU/strip) and tubes of nutrient broth with pH indicator.	20 strips 20 tubes	91072
<b>STRIP CONTROL PUM E5</b> Strips with spores of <i>B. pumilus</i> (1-5x10 <sup>5</sup> CFU/strip) and tubes of nutrient broth with pH indicator.	20 strips 20 tubes	91073 #
<b>STRIP CONTROL PUM E4</b> Strips with spores of <i>B. pumilus</i> (1-5x10 <sup>4</sup> CFU/strip) and tubes of nutrient broth with pH indicator. 20 tubes	20 strips 20 tubes	91074 #
<b>STRIP CONTROL PUM E3</b> Strips with spores of <i>B. pumilus</i> (1-5x10 <sup>3</sup> CFU/strip) and tubes of nutrient broth with pH indicator. 20 tubes	20 strips 20 tubes	91075 #