

Strips for steam sterilization

Product with minimum request



Description	Packaging	Ref.
STRIP CONTROL GST E6 Strips with spores of <i>G. stearothermophilus</i> (1-5x10 ⁶ CFU/strip) and tubes of nutrient broth with pH indicator.	20 strips 20 tubes	91055
	100 strips	91155
STRIP CONTROL GST E5 Strips with spores of <i>G. stearothermophilus</i> (1-5x10 ⁵ CFU/strip) and tubes of nutrient broth with pH indicator.	20 strips 20 tubes	91056 #
STRIP CONTROL GST E4 Strips with spores of <i>G. stearothermophilus</i> (1-5x10 ⁴ CFU/strip) and tubes of nutrient broth with pH indicator.	20 strips 20 tubes	91057 #
STRIP CONTROL GST E3 Strips with spores of <i>G. stearothermophilus</i> (1-5x10 ³ CFU/strip) and tubes of nutrient broth with pH indicator.	20 strips 20 tubes	91058 #